



SURGICAL SUPPLY CO., INC.

Once again
FASTER is not
always **BETTER**

RCM⁶
Resorbable Collagen MEMBRANE



**You CAN'T LOSE
with ACE's RCM⁶
6 Month RESORBABLE
COLLAGEN MEMBRANE**

- Fully Resorbable
- Contours Well to Defects
- Easily Cut, Shaped & Applied
- Cost Effective

www.acesurgical.com

ACE Surgical Supply Company, Inc. • 1034 Pearl Street, Brockton, MA 02301

RCM⁶ RESORBABLE COLLAGEN MEMBRANE

RESORBABLE COLLAGEN MEMBRANE is intended for use in oral surgical procedures as a resorbable material for placement to aid in wound healing in the areas of dental implants, bone defects, periodontal defects and ridge augmentation.

ACE RCM⁶™ – With its defined geometry, in vivo stability, permeability and mechanical strength, is a white, nonfriable membrane matrix engineered from highly purified Type I Collagen. ACE RCM⁶™ is packaged sterile, non-pyrogenic, and is for single use. Its resorbability eliminates the need for a second surgery.

FEATURES

- Guides healing of bone and surrounding tissue.
- Unique fiber orientation to provide high tensile strength for *in situ* stability of membrane.
- Cell occlusive, designed to retard and/or prevent epithelium down growth.

ADVANTAGES

- Longer resorption time, for sustained function in bone repair.
- Higher mechanical strength for membrane stabilization.
- Each piece can be easily trimmed during surgery for a final fit.
- May be placed with either side against the tissue or bone.

TYPICAL USES

- Ridge Augmentation
- Bone Regeneration around Implants
- Sinus Window
- Extraction Sites
- Periodontal Defects

Number	Description	Membrane size	Thickness	Price
509-6000	RCM ⁶ ™ Resorbable Collagen Mambrene	15 x 20mm	.3mm	\$89.99
509-3000	RCM ⁶ ™ Resorbable Collagen Mambrene	20 x 30mm	.3mm	109.99
509-8000	RCM ⁶ ™ Resorbable Collagen Mambrene	30 x 40mm	.3mm	164.99

ACE PGA-RAPID ABSORBING* SUTURES EQUIVALENT TO VICRYL RAPIDE®**

Number	Size	Description	Needle	Length	Quantity	Price
003-9907	3-0	Undyed Braided PGA-RA	DSM18	18"	12/box	\$79.29
003-9906	4-0	Undyed Braided PGA-RA	DSM18	18"	12/box	79.29
003-9900	4-0	Undyed Braided PGA-RA	HR17 (Taper)	27"	12/box	68.79
003-9904	5-0	Undyed Braided PGA-RA	DSM13	18"	12/box	79.29



*Maintains 50% tensile strength at 5 days; 0% at 14 days.

**Vicryl® and Vicryl Rapide® are Registered Trademarks of Ethicon, Inc.

BioMend® and BioMend® Extend are registered trademarks of Integra LifeSciences Corp. for Zimmer Dental.



Resorbable Collagen MEMBRANE



www.acesurgical.com

www.acesurgical.com